

Security Advisory

“BadAlloc” Vulnerability Assessment and Potential Product Impact Statement

Ref Doc ID	Version	Release Date	Advisory Status	Related CVE(s)	Operational Risk
079-0242-00	B	4 June 2021	ACTIVE	See ICS Advisory (ICSA-21-119-04) for details.	Low

1. VULNERABILITY

Spacelabs Healthcare has been made aware of recently published vulnerabilities collectively known as “BadAlloc” that details vulnerabilities found in multiple Real-Time Operating Systems (RTOS) and supporting libraries. Cybersecurity & Infrastructure Security Agency (CISA) issued an advisory to provide notice of the reported vulnerabilities and identify baseline mitigations for reducing risks to these and other cybersecurity attacks.

The 23 vulnerabilities identified across 25 RTOS by the security researchers from Microsoft research groups reported that successful exploitation of these vulnerabilities could result in unexpected behavior such as a crash or a remote code injection/execution.

Spacelabs has conducted an assessment to identify the potential impact on our products. Our assessment has found that Spacelabs Patient monitors including Qube, Qube Mini, Xprezzon, and UVSL Monitors are impacted by one of the identified vulnerabilities.

Spacelabs utilizes Wind River VxWorks that was identified as one of the RTOS affected by BadAlloc for the Patient Monitors. More details can be found on Wind River’s website; links are available in the reference section.

Security Impact:

Security Impact	Notes
Availability	Technical Impact: <i>DoS: Crash, Exit, or Restart; DoS: Resource Consumption (CPU); DoS: Resource Consumption (Memory); DoS: Instability</i>
Integrity	Technical Impact: <i>Modify Memory</i>
Confidentiality	Technical Impact: <i>Execute Unauthorized Code or Commands; Bypass Protection Mechanism</i>

2. OPERATIONAL RISK ASSESSMENT SUMMARY

Spacelabs considers the operational risk to its devices from a BadAlloc attack to be low. In our review there is only a limited chance that a device could be compromised due to the fact that a properly configured firewall in the hospital network would mitigate the likelihood of an attack from outside the organization and since remote code execution would require that a bad actor already be positioned within the hospital network as a result of a prior attack.

Common hospital protocol involves the use of multiple parameters and physical assessment of patients in determining treatment. For this reason, in the event that a Spacelabs monitor is compromised due to an attack exploiting BadAlloc vulnerabilities there is less likelihood of incorrect treatment. There are also typically redundant alarm systems in place to decrease the likelihood of delay in treatment.

Spacelabs is working with the vendor to develop software patches that will be available to keep an attack from using the identified operating system vulnerabilities. In the meantime, customers are encouraged to ensure their facility maintains a secure network firewall to protect its products and data.

3. RECOMMENDATIONS

Spacelabs recommends the following defenses and mitigations be applied to the enterprise environment.

- Minimize network exposure for all patient monitoring devices such as monitors (Qube, Qube mini, Xprezzon and UVSL) with the use of network segmentation, placement of these devices behind hospital firewalls, and ensuring that they are not accessible from the Internet.
- Block suspicious external IP addresses at the hospital firewalls. Monitor traffic internally for unusual behavior.
- When remote access is required, use secure methods, such as Virtual Private Networks (VPNs), recognizing VPNs may have vulnerabilities and should be updated to the most current version available. Also recognize VPN is only as secure as its connected devices.
- Implement defense-in-depth within the enterprise environment consisting of tools such as Intrusion Detection/Prevention Systems (IDS/IPS), firewalls, and network access control (NAC).
- Implement and maintain an anti-malware solution (also called “anti-virus”) and an endpoint detection and response (EDR) solution.
- Disable remote access services and protocols such as Remote Desktop Protocol (RDP) unless needed. Monitor and restrict remote access usage on a least-privilege basis.
- Have backup and restore processes and procedures in place for disaster recovery and incident response.
- Monitor and maintain account provisioning and access control based on the principle of least privilege.

4. EXAMINATION OF SPACELABS PRODUCTS

ASSESSMENT OF SPACELABS PRODUCTS

In response to the publication of these vulnerabilities, Spacelabs has conducted an assessment to identify devices potentially at risk to this set of vulnerabilities. Please note information is subject to change as the situation evolves.

Patient Monitoring and Connectivity Products

Product	Host Operating System	Impact Assessment
XprezzNet 96190	Windows Server 2012 R2, Windows Server 2016	Not impacted.
Intesys Clinical Suite (ICS)	Windows Server 2012 R2, Windows Server 2016, Windows Server 2019	Not impacted.
Xhibit Telemetry Receiver (XTR) 96280	Windows Embedded Standard 7 SP1	Not impacted.
Xhibit 96102 / XC4 96501	Windows Embedded Standard 7 SP1	Not impacted.
Bedside Monitors <ul style="list-style-type: none"> • Xprezzon 91393 • Qube 91390 • Qube Mini 91389 • Ultraview SL 91367, 91369, 91370, and 91387 	VxWorks 6.6	Impacted. There is currently no patch for the version of VxWorks running in the monitors. Spacelabs will work with RTOS developers and will issue a revision to this document when the patch is available.

Diagnostic Cardiology Products

Product	Host Operating System	Impact Assessment
Sentinel	Windows 7 & 10, Windows Server 2012 R2, Windows Server 2016, Windows Server 2019	Not impacted.
Pathfinder SL	Windows 7, Windows 10	Not impacted.
Lifescreeen Pro	Windows 10	Not impacted.
Lifecard CF	No OS	Not impacted.
EVO	No OS	Not impacted.
Eclipse Pro	No OS	Not impacted.
CardioExpress SL6A / SL12A	Embedded OS (uC/OS II V2.84)	Pending Review.
CardioExpress SL18A	Embedded OS (Linux Kernel 2.6.35.3)	Not impacted.

SPACELABS HEALTHCARE

Product	Host Operating System	Impact Assessment
ABP <ul style="list-style-type: none"> OnTrak 90217A 90207 	No OS	Not impacted.

Safe N Sound

Product	Host Operating System	Impact Assessment
Spacelabs Cloud	Varies	Not impacted.
SafeNSound	Not applicable	Not impacted.

5. Additional Resources

- Microsoft Security Response Center - <https://msrc-blog.microsoft.com/2021/04/29/badalloc-memory-allocation-vulnerabilities-could-affect-wide-range-of-iot-and-ot-devices-in-industrial-medical-and-enterprise-networks/>
- CISA ICS Advisory ICSA-21-119-04 - <https://us-cert.cisa.gov/ics/advisories/icsa-21-119-04>
- Wind River - <https://support2.windriver.com/index.php?page=cve&on=view&id=CVE-2020-28895>
- Wind River (permission required): <https://support2.windriver.com/index.php?page=defects&on=view&id=V7LIBC-1327>
- Spacelabs Patch Qualification Customer Portal - https://www.spacelabshealthcare.com/products/security/patch-test-reports-access-form/?redirect_to=%2Fproducts%2Fsecurity%2Fpatch-test-reports%2F

6. Document History

Version	Release Date	Purpose
Rev B	4 June 2021	BadAlloc Vulnerability Assessment and Potential Product Impact Statement

7. Terms of Use

The information in this document is subject to change without notice. In no event will Spacelabs or any of its suppliers be liable for direct, indirect, special, incidental, or consequential damages of any nature or kind arising from the use of this document, even if Spacelabs or its suppliers have been advised of the possibility of such damages.

This document contains confidential and proprietary language and may not be reproduced or shared with a third party without written permission from Spacelabs. All rights to registrations and trademarks reside with their respective owners.