

SUPPLIES & ACCESSORIES

NON-INVASIVE BLOOD PRESSURE CUFFS



TRULINK® NYLON REUSABLE BLOOD PRESSURE CUFFS



TruLink Nylon Reusable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Color coded for accurate sizing
- Manufactured without PVC/DEHP
- Provides a fluid barrier
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® SOFT DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Soft Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Vinyl Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Comfortable soft vinyl
- Easy to clean with PDI cloths
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1122-10	Single Tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single Tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single Tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single Tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single Tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single Tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single Tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in three adult/pediatric sizes.

- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Yellow Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1135-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1145-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1165-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® NEONATAL NON-PVC BLOOD PRESSURE CUFFS



TruLink Neonatal Non-PVC Blood Pressure Cuffs are available in five neonatal sizes for accurate measurements.

- Color coded for accurate sizing
- Genuine low profile Velcro® closures
- Manufactured without natural latex
- Fabric fold is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Neonatal Non-PVC BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



Mini Male Luer
(ML316)

715-2114-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
715-2124-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
715-2134-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
715-2144-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
715-2154-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box

Single Tube Hose



714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube, 91220	Neonatal	275 cm/9 ft	1 Each

SPACELABS NYLON REUSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Nylon Reusable Blood Pressure Cuffs are available in both regular and long sizes.

- Long size cuffs for children, small adults, adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



**Metal Bayonet
(MHP)**

714-0020-00	Single tube, MHP	Child	12 - 19 cm	1 Each
714-0070-00	Single tube, MHP	Child Long	12 - 19 cm	1 Each
714-0021-00	Single tube, MHP	Small Adult	17 - 25 cm	1 Each
714-0071-00	Single tube, MHP	Small Adult Long	17 - 25 cm	1 Each
714-0023-00	Single tube, MHP	Adult	23 - 33 cm	1 Each
714-0073-00	Single tube, MHP	Adult Long	23 - 33 cm	1 Each
714-0025-00	Single tube, MHP	Large Adult	31 - 40 cm	1 Each
714-0075-00	Single tube, MHP	Large Adult Long	31 - 40 cm	1 Each
714-0027-00	Single tube, MHP	Thigh	38 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS VINYL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Vinyl Disposable Blood Pressure Cuffs are available in seven sizes for accurate measurement.

- Long size cuffs for adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

714-1020-00	Single tube, HP	Child	12 - 19 cm	5/Pack
714-1021-00	Single tube, HP	Small Adult	17 - 25 cm	5/Pack
714-1023-00	Single tube, HP	Adult	23 - 33 cm	5/Pack
714-1073-00	Single tube, HP	Adult Long	23 - 33 cm	5/Pack
714-1025-00	Single tube, HP	Large Adult	31 - 40 cm	5/Pack
714-1075-00	Single tube, HP	Large Adult Long	31 - 40 cm	5/Pack
714-1027-00	Single tube, HP	Thigh	38 - 50 cm	5/Pack

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS SOFT CUFF NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Soft Cuff Neonatal BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



714-1028-01	Single tube, 3/16" ML	Size 1	3 - 6 cm	10/Pack
714-1029-01	Single tube, 3/16" ML	Size 2	4 - 8 cm	10/Pack
714-1030-01	Single tube, 3/16" ML	Size 3	6 - 11 cm	10/Pack
714-1032-01	Single tube, 3/16" ML	Size 4	7 - 13 cm	10/Pack
714-1033-01	Single tube, 3/16" ML	Size 5	8 - 15 cm	10/Pack



Single Tube Hose

714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube hose, 91220	Neonatal	275 cm/9 ft	1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Non-Invasive Blood Pressure Hose

Single Tube Hose

For use with Spacelabs Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

Double Tube to Single Tube Hose

For use with Command Modules 90496/91496



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------

Hose Connectors

Slip Fit Connector

Connectors fit Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.



714-0038-00	3mm ID Adult NIBP hoses			1 Each
714-0047-00	4mm ID Adult NIBP hoses			1 Each
714-0046-00	5mm ID Adult NIBP hoses			1 Each
714-0045-00	6mm ID Adult NIBP hoses			1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

NIBP Cuff Conversion Accessories

Hose Conversion Kit



714-0043-00	Dual hose conversion kit Use with TruLink and Spacelabs cuffs and dual hose blood pressure monitors incorporating thread type fittings	Adult		1 Each
-------------	---	-------	--	--------

Y-Adapters



714-0056-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single-tube 3mm ID hose and monitor/modules	Adult		5/Pack
-------------	--	-------	--	--------



714-0057-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single tube NIBP hoses and module/monitors	Adult		5/Pack
-------------	---	-------	--	--------

NIBP accessories for use with Spacelabs and TruLink branded neonatal cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Neonatal BP Cuff Accessories

Single Tube Hose

For use with Command Modules 90496/91496 and mCare Monitor 91220



714-0019-01	90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	91220	Neonatal	275 cm/9 ft	1 Each

Double Tube Cuff to Single Tube Hose

For use with Command Modules 90496/91496



714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
-------------	-------------	----------	-------------	--------

Metal Cuff Fitting

For use with NIBP hose 714-0019-01



714-0053-10	Hose fitting, metal	Neonatal		10/Pack
-------------	---------------------	----------	--	---------

Plastic Cuff Fittings

For use with NIBP Hoses 714-0019-01 & 714-0019-12



714-0052-50	Cuff fitting, plastic	Neonatal		50/Pack
714-0054-50	Cuff fitting, plastic	Infant		50/Pack

Hose Connectors

714-0057-10	NIBP non-luer hose metal connector Connects neonatal cuffs to non-Spacelabs monitors with a 3mm ID hose	Neonatal		10/Pack
-------------	--	----------	--	---------



ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Bayonet (HP)

US0814HP	Single tube, HP	Infant	8 - 14 cm	Each
US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
US1320HP	Single tube, HP	Child	13 - 20 cm	Each
US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
US1826HP	Single tube, HP	Small Adult	18 - 26 cm	Each
US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
US2635HP	Single tube, HP	Adult	26 - 35 cm	Each
US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
US2938HP	Single tube, HP	Adult Long	29 - 38 cm	Each
US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
US3242HP	Single tube, HP	Large Adult	32 - 42 cm	Each
US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
US3544HP	Single tube, HP	Large Adult Long	35 - 44 cm	Each
US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
US4250HP	Single tube, HP	Thigh	42 - 50 cm	Each
US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Metal Bayonet (MHP)

US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag



Male Luer (ML)

US0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag
US1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag
US1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag
US2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
US2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
US3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
US3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
US4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag



Screw (DM)

US0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
US2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Twist Lock (MQ)

US0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
US1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag



Female Luer (FL)

US0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
US1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
US1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
US2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
US2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
US3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
US3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
US4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag



Male Twist Lock (MQML)

US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
US1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
US2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
US4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



US0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
US1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
US1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
US4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag



US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
US1320MB	Single tube, MB w/gauge	Child	13 - 20 cm	1 Each
US1826MB	Single tube, MB w/gauge	Small Adult	18 - 26 cm	1 Each
US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
US3242MB	Single tube, MB w/gauge	Large Adult	32 - 42 cm	1 Each
US3544MB	Single tube, MB w/gauge	Large Adult Long	35 - 44 cm	1 Each
US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Double Tube



Screw (DM)

UD0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag



Twist Lock (MQ)

UD0814HP-05	Double tube, HP	Infant	8 - 14 cm	5/Bag
UD1320HP-05	Double tube, HP	Child	13 - 20 cm	5/Bag
UD1826HP-05	Double tube, HP	Small Adult	18 - 26 cm	5/Bag
UD2635HP-05	Double tube, HP	Adult	26 - 35 cm	5/Bag
UD2938HP-05	Double tube, HP	Adult Long	29 - 38 cm	5/Bag
UD3242HP-05	Double tube, HP	Large Adult	32 - 42 cm	5/Bag
UD3544HP-05	Double tube, HP	Large Adult Long	35 - 44 cm	5/Bag



Twist Lock (MQ)

UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag
UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
UD1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag
UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
UD2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag
UD3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag
UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag



Male & Female MQ (MQMF)

UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
UD1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
UD1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
UD3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
UD3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs



No Fitting

Double Tube					
UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag	
UD1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag	
UD1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag	
UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag	
UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag	
UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag	
UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag	
UD4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag	



Manual Bulb (MB)

Double Tube					
UD0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag	
UD1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag	
UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag	
UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag	
UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag	
UD3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag	
UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag	
UD4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag	

ULTRACHECK® CURVE™ REUSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Manufactured without DEHP
- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty	List Price
------------	-------------	------	-------	-----	------------

UltraCheck Curve Reusable BP Cuffs

Single Tube



BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854DM	Single tube, DM	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854FL	Single tube, FL	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854HP	Single tube, HP	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854MHP	Single tube, metal HP	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854ML	Single tube, ML	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00
BRUS3854MQ	Single tube, MQ	Extra Large Adult	38 - 54 cm	1 Each	\$ 25.00

Double Tube



BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854HP	Double tube, HP	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854ML	Double tube, ML	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854MQMF	Double tube, MQMF	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00
BRUD3854MQML	Double tube, MQML	Extra Large Adult	38 - 54 cm	1 Each	\$ 29.00

Manual Cuffs



BRUS3854MB	Single tube, manual w/ gauge	Extra Large Adult	38 - 54 cm	1 Each	\$ 52.64
BRUD3854MB	Double tube, manual	Extra Large Adult	38 - 54 cm	1 Each	\$ 38.08

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Bayonet (HP)

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box



Male Luer (ML)

ST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
ST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
ST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box
ST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
ST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box
ST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
ST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box
ST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag
ST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

ST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
ST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST1826MQ-20	Single tube, MQ	Small Adult	18 - 26 cm	20/Box
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
ST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
ST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
ST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
ST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
ST3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
ST3544MQ-20	Single tube, MQ	Large Adult Long	35 - 44 cm	20/Box
ST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
ST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Screw (DM)

ST0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
ST1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
ST1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
ST2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
ST2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
ST3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
ST3544DM-05	Single tube, DM	Large Adult Long	35 - 44 cm	5/Bag
ST4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS






SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube					
 Female Luer (FL)	ST0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
	ST1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
	ST1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
	ST2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
	ST2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
	ST3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
	ST3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
	ST4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag
 Male Twist Lock (MQML)	ST0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
	ST1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
	ST1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
	ST2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
	ST2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
	ST3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
	ST3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
	ST4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag
 No Fitting	ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
	ST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
	ST1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	ST2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	ST3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Double Tube					
DT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag	
DT0814DM-20	Double tube, DM	Infant	8 - 14 cm	20/Box	
DT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag	
DT1320DM-20	Double tube, DM	Child	13 - 20 cm	20/Box	
DT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag	
DT1826DM-20	Double tube, DM	Small Adult	18 - 26 cm	20/Box	
DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag	
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box	
DT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag	
DT2938DM-20	Double tube, DM	Adult Long	29 - 38 cm	20/Box	
DT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag	
DT3242DM-20	Double tube, DM	Large Adult	32 - 42 cm	20/Box	
DT3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag	
DT3544DM-20	Double tube, DM	Large Adult Long	35 - 44 cm	20/Box	
DT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag	
DT4250DM-20	Double tube, DM	Thigh	42 - 50 cm	20/Box	



Screw (DM)

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Double Tube					
DT0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag	
DT0814MQ-20	Double tube, MQ	Infant	8 - 14 cm	20/Box	
DT1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag	
DT1320MQ-20	Double tube, MQ	Child	13 - 20 cm	20/Box	
DT1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag	
DT1826MQ-20	Double tube, MQ	Small Adult	18 - 26 cm	20/Box	
DT2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag	
DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box	
DT2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag	
DT2938MQ-20	Double tube, MQ	Adult Long	29 - 38 cm	20/Box	
DT3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag	
DT3242MQ-20	Double tube, MQ	Large Adult	32 - 42 cm	20/Box	
DT3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag	
DT3544MQ-20	Double tube, MQ	Large Adult Long	35 - 44 cm	20/Box	
DT4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag	
DT4250MQ-20	Double tube, MQ	Thigh	42 - 50 cm	20/Box	



Twist Lock (MQ)

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Double Tube



Male & Female MQ (MQMF)

DT0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
DT1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
DT1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
DT2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
DT2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
DT3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
DT3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
DT4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag



No Fitting

DT0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
DT1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag
DT1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
DT2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
DT3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
DT3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
DT4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag



Manual Bulb (MB)

DT0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
DT1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag
DT1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
DT2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
DT3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
DT3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
DT4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® CURVE™ DISPOSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty	List Price
------------	-------------	------	-------	-----	------------

SoftCheck Curve Disposable BP Cuffs

Single Tube						
BRST3854HP-05	Single tube, HP	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854DM-05	Single tube, DM	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854FL-05	Single tube, FL	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854ML-05	Single tube, ML	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854MQ-05	Single tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854MQML-05	Single tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag	\$ 50.00	
Double Tube						
BRDT3854HP-05	Double tube, HP	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854DM-05	Double tube, DM	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854FL-05	Double tube, FL	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854ML-05	Double tube, ML	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854MQ-05	Double tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854MQMF-05	Double tube, MQMF	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854MQML-05	Double tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
BRDT3854-05	Double tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag	\$ 60.00	
Manual Cuff						
BRDT3854MB-05	Double tube, MB	Extra Large Adult	38 - 54 cm	5/Bag	\$ 91.43	



SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs



Bayonet (HP)

Single Tube				
VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
VST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
VST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
VST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
VST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
VST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
VST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
VST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube



Male Luer (ML)

VST0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag
VST0814ML-20	Single tube, ML	Infant	8 - 14 cm	20/Box
VST1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag
VST1320ML-20	Single tube, ML	Child	13 - 20 cm	20/Box
VST1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag
VST1826ML-20	Single tube, ML	Small Adult	18 - 26 cm	20/Box
VST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
VST2635ML-20	Single tube, ML	Adult	26 - 35 cm	20/Box
VST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	20/Box
VST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box
VST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
VST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box
VST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
VST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box
VST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag
VST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

VST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
VST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
VST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
VST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
VST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
VST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
VST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
VST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
VST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Male Twist Lock (MQML)

VST0814MQML-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQML-20	Single tube, MQML	Child	13 - 20 cm	20/Box
VST1826MQML-20	Single tube, MQML	Small Adult	18 - 26 cm	20/Box
VST2635MQML-20	Single tube, MQML	Adult	26 - 35 cm	20/Box
VST2938MQML-20	Single tube, MQML	Adult Long	29 - 38 cm	20/Box
VST3242MQML-20	Single tube, MQML	Large Adult	32 - 42 cm	20/Box
VST3544MQML-20	Single tube, MQML	Thigh	35 - 44 cm	20/Box
VST4250MQML-20	Single tube, MQML	Thigh	42 - 50 cm	20/Box



No Fitting

VST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
VST1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
VST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
VST3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
VST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS





SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Double Tube					
 Screw (DM)	VDT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
	VDT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	VDT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	VDT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	VDT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	VDT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	VDT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
 No Fitting	VDT0814-20	Double tube, no fitting	Infant	8 - 14 cm	20/Box
	VDT1320-20	Double tube, no fitting	Child	13 - 20 cm	20/Box
	VDT1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	VDT2635-20	Double tube, no fitting	Adult	26 - 35 cm	20/Box
	VDT2938-20	Double tube, no fitting	Adult Long	29 - 38 cm	20/Box
	VDT3242-20	Double tube, no fitting	Large Adult	32 - 42 cm	20/Box
	VDT4250-20	Double tube, no fitting	Thigh	42 - 50 cm	20/Box

SOFTCHECK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Yellow Vinyl Disposable BP Cuffs



Bayonet (HP)

Single Tube					
YVST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag	
YVST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag	
YVST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag	

SOFTCHECK® FABRIC NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Neonatal Disposable BP Cuffs

Single Tube



Mini Male Luer
(ML316)

NN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
NN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
NN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
NN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
NN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box



Male Luer (ML)

NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Bayonet (HP)

NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
NN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
NN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



Male Luer (ML)

NN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
NN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
NN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
NN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS








SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs

Single Tube					
 ML316	VNN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box
 SnapQuik (SQ)	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 Male Luer (ML)	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Bayonet (HP)	VNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
	VNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
	VNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
	VNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
	VNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box
 Twist Lock (MQ)	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.

- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs



Male Luer (ML)

Double Tube					
VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box	
VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box	
VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box	
VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box	
VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box	

SOFTCHECK FREE® NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Free Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Created with extremely soft material
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Free Non-PVC Neonatal Disposable BP Cuffs

Single tube



PNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
PNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
PNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
PNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
PNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box



PNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
PNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
PNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
PNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
PNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



PNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
PNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
PNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
PNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
PNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



PNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
PNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
PNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
PNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
PNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
------------	-------------	-----

Blood Pressure Cuff Fittings & Adapters

Used in conversions to HP Cuffs

ZFT-HP-20	Male Bayonet Connector for Disposable Cuffs	20/Pack
ZFT-HP532-20	Male Bayonet Connector for Reusable Cuffs	20/Pack
ZFT-QC8030	Male Metal Bayonet Connector, 5/32" barb	1 Each
ZFT-QC8220	Female Metal Bayonet Connector, 1/8" barb	1 Each
ZFT-QC8230	Female Metal Bayonet Connector, 5/32" barb	1 Each
ZFT-BPF230	Female Plastic Bayonet Connector, 1/8" barb	1 Each
ZFT-BPF240	Female Plastic Bayonet Connector, 5/32" barb	1 Each
ZHB-BPCOILADPRX1	Coiled Tube Manual Bulb Adapter w/metal female bayonet fitting	1 Each
ADP-182YBPF	Y Adapter w/1/8" barb fittings & plastic female bayonet connector	1 Each
ADP-182YRX4	Y Adapter w/1/8" barb fittings & metal female bayonet connector	1 Each
ADP-5322YBPF	Y Adapter w/5/32" barb fittings & plastic female bayonet connector	1 Each
ADP-5322YRX4	Y Adapter w/5/32" barb fittings & metal female bayonet connector	1 Each
ADP-DM2YRX4	Y Adapter, Metal Female HP to 2 DM	1 Each
ZFT-HPDMAD-05	Y Adapter /Male Bayonet to 2 Male Screw	5/Pack
ZFT-ML2YFL	Y Adapter w/Female Luer to 2 Male Luers	1 Each

Used in conversions to HP Cuffs-Straight

ZFT-RX1DM	5" Tube w/Metal Female Bayonet and Screw Fitting	1 Each
ZFT-RX1FL	5" Tube w/Metal Female Bayonet and Female Luer	1 Each
ZFT-RX1ML	5" Tube w/Metal Female Bayonet and Male Luer	1 Each
ZFT-RX1MQ	5" Tube w/Metal Female Bayonet and Twist-Lock Fitting	1 Each
ZFT-RX1MQML	5" Tube w/Metal Female Bayonet and Male Twist-Lock Fitting	1 Each

Used in conversions to HP Cuffs-Ys

ZFT-DM2YRX1	Y Adapter, Metal Female Bayonet to 2 Screw fittings	1 Each
ZFT-MQ2YRX1	Y Adapter w/Metal Female Bayonet to 2 Twist Lock fittings	1 Each
ZFT-MQML2YRX1	Y Adapter w/Metal Female Bayonet to 2 Male Twist Lock fittings	1 Each
ZFT-MQMF2YRX1	Y Adapter w/Metal Female Bayonet to one female Twist Lock and one male Twist Lock fitting	1 Each

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
Blood Pressure Cuff Fittings & Adapters		
Miscellaneous		
ZFT-CS901	Brass Male Screw Fitting	1 Each
ZFT-DM-20	Female Screw Connector for Disposable Cuffs	20/Pack
ZFT-DM532-20	Female Screw Connector for Reusable Cuffs	20/Pack
ZFT-DMML-20	Nylon Male Screw Fitting	20/Pack
ZFT-FL532-20	Female Luer for Disposable Cuffs	20/Pack
ZFT-FLL240-20	Female Luer for Reusable Cuffs	20/Pack
ZFT-ML532-20	Male Slip Luer for Disposable and Reusable Cuffs	20/Pack
ZFT-MLL532-20	Male Locking Luer	20/Pack
ZFT-MQ-20	Female Twist-Lock Fitting for Disposable and Reusable Cuffs	20/Pack
ZFT-MQMF2YHP	Y Adapter - Twist Lock and Male Twist Lock to Male Bayonet	1 Each
ZFT-MQML-20	Male Twist-Lock Fitting for Disposable Cuffs	20/Pack
ZFT-MQML5-20	Male Twist-Lock Fitting for Reusable Cuffs	20/Pack
ZGA-BPGBHB	BP Gauge w/Black 2 oz. Handbulb	1 Each
ZGA-BPHBRX1	Aneroid Gauge BPC Inflator, 8 in. tubing w/Metal Female Bayonet	1 Each
ZHB-4FBPADPRX1	4' single tube Manual Bulb BP Cuff Adapter w/Metal Female Bayonet	1 Each
ZHB-HBTHM	Black 2oz Handbulb w/Metal Thumbwheel	1 Each