



*Xhibit® XC4 (96501)*

The Spacelabs Healthcare Xhibit® XC4 (96501) provides central monitoring of patient data in the patient vicinity on a medical grade all-in-one computer with display. You can monitor and manage alarms on one to four patients at any location with connectivity to the Spacelabs Healthcare patient monitor network. This four-bed central station at the point-of-care provides clinicians with the ability to monitor other beds, including telemetry patients. The Xhibit XC4 provides full interactive capability to admit, discharge, and modify patient information and to update parameter settings.

## Features

Patient Capacity	<ul style="list-style-type: none"> <li>Shows all parameters monitored for up to four patients on a high resolution 15.6-inch touchscreen display</li> <li>Full Bedside View shows all parameters for any patient, while maintaining visibility for all monitored patients</li> </ul>
User Interface	Intuitive controls on the touchscreen interface provide simplified patient management. Keyboard and mouse are additional options for use.
Printing	Strips and reports print directly from the Xhibit XC4. Receives and prints strips from both bedside monitors and telemetry transmitters.
Trend Review	View any parameter trend data for up to 96 hours in graphic or tabular form

## Product Specifications

CPU	Intel Core i5-4300U Processor
Memory	8 GB total
Storage	256 GB Solid State Drive (SSD)
Display	1920 x 1080 native resolution configured to 1680 x 1050
Size	39.6 cm (15.6 inch) widescreen with Thin-Film Transistor (TFT) LCD
Operating System	Windows 10 Enterprise LTSC
Mains Power	<ul style="list-style-type: none"> <li>100 to 240 VAC</li> <li>1.2 to 0.5 A</li> <li>47 to 63 Hz</li> </ul>
Touchscreen	<ul style="list-style-type: none"> <li>Capacitive Touchscreen (P-cap) that supports gloved finger touch</li> <li>Surface P-cap is anti-reflection (AR) and anti-glare (AG)</li> </ul>

## Environmental Requirements

Storage	
Temperature	-20° to 65° C (-4° to 167° F)
Humidity	5 to 95% (noncondensing)
Altitude	up to 12,192 meters (40,000 feet)
Operating	
Temperature	0° to 40° C (32° to 104° F)
Humidity	10 to 95% at 40° C (noncondensing)
Altitude	up to 3,000 meters (10,000 feet)

## Physical Characteristics

### Dimensions

Height	271 mm (10.67 inches)
Width	405 mm (15.94 inches)
Depth	59.8 mm (2.35 inches)
Weight	4.65 kg (10.25 lbs)

### Accessories

DC Power Supply	P/N 119-0596-00
Printers	The printers are office equipment and not for use in the patient vicinity. <ul style="list-style-type: none"> <li>• HP LaserJet M507, 110V (P/N 010-1971-02)</li> <li>• HP LaserJet M507, 220V (P/N 010-1886-02)</li> </ul>

## Classification

MDD	Class IIb medical device, rule 10 <ul style="list-style-type: none"> <li>• EN60601-1</li> <li>• EN60601-1-6</li> <li>• EN60601-1-8</li> </ul>
-----	---

## Licensing

Patient Zone License	P/N 910-6501-01
----------------------	-----------------

## Compatibility

The Xhibit XC4 is compatible with Spacelabs Healthcare Ultraview SL<sup>®</sup>, Qube<sup>®</sup>, Qube<sup>®</sup> Mini, and Xprezzon<sup>®</sup> bedside monitors, as well as the Xhibit<sup>®</sup> Telemetry Receiver. The Xhibit XC4 is also compatible with selected Flexports.

*Note:*

*Not all products are available in all locales. Check with your local Spacelabs Healthcare representative.*

## Documentation

This product ships with a complete set of comprehensive documentation.

For a full list of available supplies and accessories, refer to the *Spacelabs Healthcare Supplies and Accessories Catalog* at <https://www.spacelabshealthcare.com/products/supplies>.

## Regulatory Approvals



UL certified. Meets IEC 60601-1; ANSI/AAMI ES 60601-1; and CSA C22.2 No. 60601-1 for use in the patient vicinity.



CE marked in accordance with the Medical Device Directive 93/42/EEC.

*We are continuously improving our products. Specifications are subject to change without notice. Product images are provided for general reference. This product may not be available for sale in all countries.*

*Please refer to <https://www.spacelabshealthcare.com/trademarks> for a full listing of Spacelabs Healthcare trademarks. Other brands and product names used herein are trademarks of their respective owners.*