

2024-8

English

Spacelabs Healthcare produces clinically beneficial, cost-effective monitoring and clinical information systems. As part of our commitment to the continual development of our product line, we have released the latest version of the product listed below. This bulletin is designed to inform you of the new features or enhancements. Thank you for your continued support.

Update specifications

Updated System Components	SafeNSound Mobile 96180
Software Version	5.0.7

New features or enhancements

Mobile updates will be distributed by the hospital organization's preferred choice of MDM Solution. This should be managed by the organization.

SafeNSound Mobile SW v5.0.7 contains minor improvements mentioned below.

Other changes and updates

- Recent Number - The most recent outgoing phone number is cleared from the dial pad and stored in call history.
- Admitting on Android - The application error when attempting to scan a patient into SafeNSound using Android devices is now resolved.
- Alarms and Tasks - Alarms and Tasks header was removed from the Dashboard. This feature will return in SafeNSound v5.1.0.

SafeNSound Mobile 5.x Disabled Features

SafeNSound Mobile 5.0.6 does not currently include the following features that were present in SafeNSound Mobile version 4.x:

- Detailed Patient Data - The ability to select a patient and view corresponding patient data is disabled.
- Messaging - This feature is disabled.
- Profile - This feature is disabled.
- Support - This feature is disabled.

Supported versions

- Web client version 8.0+
- Smart Alarm Interface (SAI-96181) version 1.2+
- WebApp - 5.9.7+